

## LAB WORK 12.

**Subject: Genetic recombination. Conjugation in Escherichia coli K-12.**Session

**Purpose:** To implement Conjugation in Escherichia coli K-12.

**Objectives:**

1. Learn of how to make a Suspension hybridization.
2. Learn of how to make a Spot hybridization.
3. Perform exercises 1, 2.

Conjugation requires cell-to-cell contact for DNA to be transferred from a donor to a recipient. Bacterial conjugation is plasmid-mediated gene transfer. A plasmid that can mediate gene transfer is termed the F (fertility) plasmid. A bacterial cell containing the F plasmid is called an F<sup>+</sup> cell. A bacterial cell not containing a F plasmid is called the F<sup>-</sup> cell. A bacterial cell containing a F plasmid integrated into the bacterial DNA is termed a Hfr (high frequency of recombination) cell.

### Lab Exercise 1. Suspension hybridization.

**Methodical instructions:** Use the Met His Val Ser Strains of Escherichia coli K-12.

Strain	Sex	genotype
1	F <sup>-</sup>	Met - His <sup>-</sup> Val <sup>+</sup> Ser <sup>+</sup>
2	F <sup>+</sup>	Met + His <sup>+</sup> Val <sup>-</sup> Ser <sup>-</sup>
3	Hfr	Met + His <sup>+</sup> Val <sup>-</sup> Ser <sup>-</sup>

Make crossing between strains: F<sup>+</sup> x F<sup>-</sup> and Hfr x F<sup>-</sup>. Conjugation is carried out in liquid medium. The conjugation mix for this purpose prepares.

**Procedure:**

1. To culture of a strain 1 add culture of a strain 2.
2. To culture of a strain 1 add culture of a strain 3.
3. Place test tubes in the thermostat, 37<sup>0</sup>, 1 hour.
4. Divide Petri's cup with the minimum medium into 2 sectors.
5. Sign each.
6. From the first test tube place 100 mcl of a conjunction mix in sector 1.
7. Distribute a mix with a spatula on surface of nutrient medium.
8. From the second test tube place 100 mcl of a conjunction mix in sector 1.
9. Distribute a mix with a spatula on surface of nutrient medium.
10. Place cups in the thermostat, 37<sup>0</sup>, 48 hours.

### Lab Exercise 2. Spot hybridization.

**Methodical instructions:** Make crossing between strains: F<sup>+</sup> x F<sup>-</sup> and Hfr x F<sup>-</sup>. Transfer of a factor occurs not only by crossing in a liquid broth, but also on a cup surface with agar.

**Procedure:**

1. Divide a Petri dish with the minimum nutrient medium into 5 sectors.
2. Sign each.
3. Drop of suspension on sector 1 of a strain 1.
4. Drop of suspension on sector 2 of a strain 2.
5. Drop of suspension on sector 3 of a strain 3.
6. On sector 4, firstly a drop of a strain 1 and when it dry a drop of suspension a strain 2.

7. On sector 4, firstly a drop of a strain 1 and when it dry a drop of suspension a strain 3.
8. Petri dish put in the thermostat on 37<sup>o</sup>, 48 hours.

Write the results of experiment in the table. Make a conclusion.

#### Results of conjugation in E. Coli

Method of hybridization	Type of hybridization	Presence of recombinants
Spot	F- x F-	
	F- x Hfr	
Suspension	F- x F+	
	F- x Hfr	

#### **Equipment:**

- Broth Culture of F-, F+, Hfr strains Escherichia coli K-12
- Petri dishes with minimum nutrient medium
- Disinfectant tray
- Inoculation loop
- Burnerflame
- Spatula